



Special Tube Instruction List (Page 1)

1	2	3	4	5	6
Container Code: GRNLI	Container Code: NAGRN	Container Code: NAGRN10	Container Code: GOLD	Container Code: YELB Solution B	Container Code: YELA Solution A
Volume: 4.0 mL	Volume: 2.0 mL	Volume: 10 mL	Volume: 5.0 mL	Volume: 6.0 mL	Volume: 8.5 mL
<u>Additive:</u> Lithium Heparin	<u>Additive:</u> Sodium Heparin	<u>Additive:</u> Sodium Heparin	<u>Additive:</u> Clot activator & gel for serum separation	<u>Additive:</u> Sodium polyanethol sulfonate. Acid citrate dextrose	<u>Additive:</u> Trisodium citrate, citric acid, dextrose.
Tube Color: Green	Tube Color: Green with a semi clear cap	Tube Color: Green	Tube Color: Gold	Tube Color: Yellow	Tube Color: Yellow



Special Tube Instruction List (Page 2)

7	8	9	10	11	12
Container Code: DBLUEDTA	Container Code: DBLU	Container Code: FDP	Container Code: LPST	Container Code: GRAY	Container Code: BLU
Volume: 6 mL	Volume: 6 mL	Volume: 2 mL	Volume: 4 mL	Volume: 4 mL	Volume: 1.8 mL
<u>Additive:</u> K2 EDTA	<u>Additive:</u> Clot Activator	<u>Additive:</u> FDP Tube	<u>Additive:</u> K2 EDTA Clot activator & gel for serum separation	<u>Additive:</u> Potassium oxalate/ sodium fluoride	<u>Additive:</u> Na Citrate (For low volume blood)
Tube Color: Royal Blue	Tube Color: Royal Blue	Tube Color: Royal Blue with a yellow label	Tube Color: Lavender with a yellow ring inside	Tube Color: Gray	Tube Color: Light Blue with clear cover

