

Request for Supplies

Yankton Service Center
 501 Summit Street, Yankton, SD 57078
 Phone: 605-668-8169 OR toll free 866-256-2744
 Fax: **605-668-8168**

Client: _____ Date: _____
 Address: _____
 Person Requesting: _____ Phone: _____

SUPPLY REQUESTS MUST BE SUBMITTED IN WRITING – BY FAX, COURIER DELIVERY, ETC.

Item (unit)	QTY ORDERED	QTY SHIPPED
SPECIMEN COLLECTION:		
Vacutainer-purple K2 EDTA 3 ml (100/flat or ea)		
Vacutainer-blue Citrate NA 1.8 ml (ea)		
Vacutainer-blue Citrate NA 2.7 ml (100/flat or ea)		
Vacutainer-green PST Gel LI Heparin 4 ml (100/flat or ea)		
Vacutainer-red top 4 ml (100/flat or ea)		
Sed rate tube (ea)		
Vacutainer-gold SST Gel serum 3.5 ml (100/flat or ea)		
Vacutainer-red/black SST Gel serum 7.5 ml (100/flat or ea)		
Vacutainer-red NO Gel serum 10 ml (100/flat or ea)		
Vacutainer-dark green NA Hep 6 ml (ea)		
Vacutainer-FDP 2 ml (ea)		
Vacutainer-gray top Na Floride 6 ml (ea)		
Vacutainer-royal blue K2 EDTA 6 ml (ea)		
Vacutainer-royal blue NO additive 6 ml (ea)		
MICROBIOLOGY/VIROLOGY COLLECTION:		
Blood Culture-Aerobic green (25/flat or ea)		
Blood Culture-Anaerobic purple (25/flat or ea)		
Blood Culture-Peds yellow (ea)		
Chlam/GC DNA Amp Probe Collect Kit (female/male) (ea)		
Chlam/GC DNA Amp Probe Collect Kit - urine (ea)		
Culture swab-aerobic (white) (ea)		
Culture swab-anaerobic (blue) (ea)		
Culture swab-mini tip (green) (ea)		
Media-MAC (ea)		
Media-BAP TSA 5% SB (ea)		
Media-Choc II (ea)		
Media-Thio (ea)		
Media-MTM (ea)		
M4 Media (ea)		
Ova & Parasite 10% Formalin/PVA (ea)		
Stool Culture-Cairy Blair (C&S) (ea)		
Vaginitis DNA Probe Affirm VPIII (collect kit) (ea)		
Bio Bag (type C)		

Item (unit)	QTY ORDERED	QTY SHIPPED
MISCELLANEOUS:		
ALN Service Guide (ea)		
Requisition forms (ea)		
Fee Schedule (ea)		
Fetal Fibronectin (Spec Coll) Kits (8/bx)		
Glucola 100 g (ea)		
Printer paper		
Printer ribbon		
Supply order form		

SPECIMEN TRANSPORT:		
Biohazard Specimen Transport Bags (100/bundle)		
Transport Vials-Plastic (Amber) (ea)		
Transport Vials-Plastic (ea)		
7 ml transfer pipets (500/box)		

URINE COLLECTION:		
24-hr urine container (ea)		
Urine Culture Transport (grey top) (ea)		
UA Tube Transport (Yellow/Red Top) (ea)		
Urine Containers-Sterile Cup (ea)		

CLIENT –or - ASHH COMMENTS/OTHER:

LAB DOCUMENTATION ONLY
Order Filled By: _____
Date Filled: _____
Sent By: In-town courier Med Express Same Day Other _____